

All variables defined in `cry_ini.py`	Found in `Mtd.pref` ?	Used Where?
RawDir	True	<pre>def focus_all(...) sampleSumLists = cry_utils.get_sample_list(EXPR_FILE.basefile, samplelistTexte, EXPR_FILE.RawDir)</pre>
VanDir	True	-
VEmptyDir	True	-
SEmptyDir	True	-
CorrVanDir	True	<pre>def load_sac_eff(EXPR_FILE, NoSAC=False, Eff=True): newCalFile = EXPR_FILE.CorrVanDir + '/' + EXPR_FILE.GrpFile</pre>
GrpDir	True	-
VanFile	False	<pre>def load_sac_eff(EXPR_FILE, NoSAC=False, Eff=True): (dum, uampstotal) = cry_sample.get_data_sum(EXPR_FILE.VanFile, "Vanadium", EXPR_FILE)  def create_vana(EXPR_FILE, NoAbs=False): (dum, uampstotal) = cry_sample.get_data_sum(EXPR_FILE.VanFile, "Vanadium", EXPR_FILE)</pre>
VEmptyFile	False	<pre>def create_vana(EXPR_FILE, NoAbs=False): (dum, uampstotal) = cry_sample.get_data_sum(EXPR_FILE.VEmptyFile, "Empty", EXPR_FILE)</pre>
SEmptyFile	False	<pre>if EXPR_FILE.SEmptyFile[0] != "none": Minus(LHSWorkspace="sample", RHSWorkspace="Sempty", OutputWorkspace="sample") mtd.remove("Sempty")</pre>
AutoVan	False	-
CorrVanFile	True	<pre>def focus_all(...)  vanfil = EXPR_FILE.CorrVanFile + "-" + str(spec) + ".nxs"  def create_vana(EXPR_FILE, NoAbs=False): SaveNexusProcessed(Filename=EXPR_FILE.CorrVanFile + "_unstripped.nxs", InputWorkspace="Vanadium") SaveFocusedXYE(Filename=EXPR_FILE.CorrVanFile + "_unstripped.dat", InputWorkspace="Vanadium", SplitFiles=True)  def save_vana(EXPR_FILE): vanfil = EXPR_FILE.CorrVanFile + "-" + str(spec) + ".nxs"  def remove_bins(EXPR_FILE): SaveFocusedXYE(Filename=EXPR_FILE.CorrVanFile + "-" + str(spec) + "_unstripped.dat", InputWorkspace="Vanadium-" + str(i),</pre>

		<pre> SplitFiles=False)  def van_strip(EXPR_FILE): SaveFocusedXYE(FileName=EXPR_FILE.CorrVanFile + "-" + str(spec) + ".dat", InputWorkspace="Vanadium-" + str(i), SplitFiles=False)  def van_spline(EXPR_FILE): SaveFocusedXYE(FileName=EXPR_FILE.CorrVanFile + "-" + str(spec) + ".dat", InputWorkspace="Vanadium-" + str(i), SplitFiles=False)  def van_spline_only(EXPR_FILE): SaveFocusedXYE(FileName=EXPR_FILE.CorrVanFile + "-" + str(spec) + ".dat", InputWorkspace="Vanadium-" + str(i), SplitFiles=False) </pre>
SacEffFile		-
ExistV	False	<pre> def focus_all(...)  if EXPR_FILE.ExistV == "load":  def create_vana(EXPR_FILE, NoAbs=False): if EXPR_FILE.ExistV == 'no' and EXPR_FILE.VGrpfocus == 'van': EXPR_FILE.write_preline("ExistingV", "yes") EXPR_FILE.ExistV = "yes" </pre>
VGrpfocus	False	<pre> def focus_one(...) if EXPR_FILE.VGrpfocus == "sam" and isfirst: cry_vana.create_vana(EXPR_FILE, NoAbs)  def create_vana(EXPR_FILE, NoAbs=False): if EXPR_FILE.ExistV == 'no' and EXPR_FILE.VGrpfocus == 'van': EXPR_FILE.write_preline("ExistingV", "yes") EXPR_FILE.ExistV = "yes" </pre>
Path2DatGrpFile	False	<pre> def focus_one(...) EXPR_FILE.Path2DatGrpFile = newCalFile ... DiffractionFocussing(InputWorkspace="sample", OutputWorkspace="sample", GroupingFileName=EXPR_FILE.Path2DatGrpFile, PreserveEvents=False)  def create_vana(EXPR_FILE, NoAbs=False): if EXPR_FILE.VGrpfocus == "sam": GrpFile = EXPR_FILE.Path2DatGrpFile </pre>
VHeight	True	<pre> cry_utils.correct_abs(InputWkspc="Vanadium_corr", outputWkspc="Transmission", \ TheCylinderSampleHeight=EXPR_FILE.VHeight, \ TheCylinderSampleRadius=EXPR_FILE.VRadius, \ TheAttenuationXSection=EXPR_FILE.VAttenuationXSection, \ TheScatteringXSection=EXPR_FILE.VScatteringXSection, \ TheSampleNumberDensity=EXPR_FILE.VanaNumberDensity, \ TheNumberOfSlices=EXPR_FILE.VNumberOfSlices, \ </pre>
VRadius	True	
VAttenuationXSection	True	
VScatteringXSection	True	
VanaNumberDensity	True	
VNumberOfSlices	True	
VNumberOfAnnuli	True	
VNumberOfWavelengt	True	

hPoints		<code>TheNumberOfAnnuli=EXPR_FILE.VNumberOfAnnuli, \</code> <code>TheNumberOfWavelengthPoints=EXPR_FILE.VNumberOfWavelengthPoints, \TheExpMethod=EXPR_FILE.VExpMethod)</code>
VExpMethod	True	
VanSmooth	True	<code>def remove_bins(EXPR_FILE):</code> <code>OutputWorkspace="Vanadium-" + str(i),</code> <code>    NPoints=int(EXPR_FILE.VanSmooth))</code>
VanSplineCoef	True -> Differs	<code>def van_spline(EXPR_FILE):</code> <code>SplineBackground(InputWorkspace="Vanadium-" + str(i),</code> <code>OutputWorkspace="Vanadium-" + str(i), WorkspaceIndex=0,</code> <code>    NCoeff=int(EXPR_FILE.VanSplineCoef))</code>  <code>def van_spline_only(EXPR_FILE):</code> <code>SplineBackground(InputWorkspace="Vanadium-" + str(i),</code> <code>OutputWorkspace="Vanadium-" + str(i), WorkspaceIndex=0,</code> <code>    NCoeff=int(EXPR_FILE.VanSplineCoef))</code>
VanPeakRemove	True -> Differs	<code>def strip_the_vana(EXPR_FILE, LoadUnstrip=""):</code> <code>    if EXPR_FILE.VanPeakRemove == "interpol":</code> <code>        print " =&gt; Van Bragg-peak stripping"</code> <code>        print "Smoth Vana data with " +</code> <code>EXPR_FILE.VanSmooth + " points"</code> <code>        remove_bins(EXPR_FILE)</code> <code>    elif EXPR_FILE.VanPeakRemove == "strip":</code> <code>        print " =&gt; Van Bragg-peak stripping"</code> <code>        van_strip(EXPR_FILE)</code> <code>    elif EXPR_FILE.VanPeakRemove == "spline":</code> <code>        van_spline(EXPR_FILE)</code> <code>    elif EXPR_FILE.VanPeakRemove == "splineonly":</code> <code>        van_spline_only(EXPR_FILE)</code> <code>    else:</code> <code>        return</code> <code>save_vana(EXPR_FILE)</code>
VanPeakFile	True	-
VanPeakList	False	<code>def remove_bins(EXPR_FILE):</code> <code>for peak in EXPR_FILE.VanPeakList[spec]:</code> <code>    RemoveBins(InputWorkspace="Vanadium-" + str(i),</code> <code>OutputWorkspace="Vanadium-" + str(i), XMin=peak[0],</code> <code>        XMax=peak[1], RangeUnit="AsInput",</code> <code>Interpolation="Linear", WorkspaceIndex=0)</code>  <code>def van_spline(EXPR_FILE):</code> <code>for peak in EXPR_FILE.VanPeakList[spec]:</code> <code>    MaskBins(InputWorkspace="Vanadium-" + str(i),</code> <code>OutputWorkspace="Vanadium-" + str(i), XMin=peak[0],</code> <code>        XMax=peak[1])</code>
VanPeakWdt	False	<code>def van_strip(EXPR_FILE):</code> <code>if EXPR_FILE.VanPeakWdt[spec] != 0:</code> <code>EXPR_FILE.VanPeakWdt[spec]) + " Tol=" +</code> <code>str(EXPR_FILE.VanPeakTol[spec])</code>  <code>StripPeaks(InputWorkspace="Vanadium-" + str(i),</code> <code>OutputWorkspace="Vanadium-" + str(i),</code> <code>    FWHM=EXPR_FILE.VanPeakWdt[spec],</code> <code>Tolerance=EXPR_FILE.VanPeakTol[spec], WorkspaceIndex=0)</code> <code>SaveFocusedXYE(...)</code>
VanPeakTol	False	
CorrectSampleAbs	True	<code>def focus_one(...)</code>

SampleAbsCorrected	False	<pre> if EXPR_FILE.CorrectSampleAbs == "yes":     if EXPR_FILE.SampleAbsCorrected == False:         cry_utils.correct_abs(InputWkspc="sample",                                outputWkspc="SampleTrans", \                                TheCylinderSampleHeight=EXPR_FILE.SampleHeight, \                                TheCylinderSampleRadius=EXPR_FILE.SampleRadius, \                                TheAttenuationXSection=EXPR_FILE.SampleAttenuationXSection, \                                TheScatteringXSection=EXPR_FILE.SampleScatteringXSection, \                                TheSampleNumberDensity=EXPR_FILE.SampleNumberDensity, \                                TheNumberOfSlices=EXPR_FILE.SampleNumberOfSlices, \                                TheNumberOfAnnuli=EXPR_FILE.SampleNumberOfAnnuli, \                                TheNumberOfWavelengthPoints=EXPR_FILE.SampleNumberOfWavelengthPoints, \                                TheExpMethod=EXPR_FILE.SampleExpMethod)          EXPR_FILE.SampleAbsCorrected = True else:     ConvertUnits(...)     Divide(...)     ConvertUnits(...)     DiffractionFocussing(...) </pre>
SampleHeight	True	
SampleRadius	True	
SampleAttenuationXSection	True	
SampleScatteringXSection	True	
SampleNumberDensity	True	
SampleNumberOfSlices	True	
SampleNumberOfAnnuli	True	
SampleNumberOfWavelengthPoints	True	
SampleExpMethod	True	
LowerLambda	True	<pre> def load_sac_eff(EXPR_FILE, NoSAC=False, Eff=True):     Integration(InputWorkspace="Eff",                 OutputWorkspace="Eff", \                 RangeLower=EXPR_FILE.LowerLambda,                 RangeUpper=EXPR_FILE.UpperLambda) </pre>
UpperLambda	True -> Differs	
OutSuf	False	<pre> def focus_one(...) if EXPR_FILE.OutSuf == "":     OutputFile = join(EXPR_FILE.user, outname) else:     OutputFile = join(EXPR_FILE.user, outname + "_" + EXPR_FILE.OutSuf) # Gss rearrang4gss(OutputFile, EXPR_FILE) # Nexus rearrange 4nex(OutputFile, EXPR_FILE) </pre>
XYEDspc	False	-
XYEtof	False	-
Drange	False	<pre> def focus_one(...) cry_load.bin_bank("ResultD", EXPR_FILE.bankList, EXPR_FILE.Drange)  def bin_bank(InputArea, bankList, Drange):     for i in bankList:         Rebin(InputWorkspace=InputArea + "-" + str(i),               OutputWorkspace=InputArea + "-" + str(i),               Params=Drange[i - 1])  def sets_drange(wkspc, EXPR_FILE):     Drange = xbegin + "," + EXPR_FILE.Bining[i] + "," + xend     EXPR_FILE.Drange.append(Drange) </pre>

## Parameters:

All variables read from `mtd.pref`	Equals `=` in `mtd.pref`	Used Where
OffFile / Offsets	hrpd_new_072_01_corr.cal	-
GrpFile / Grouping	hrpd_new_072_01_corr.cal	<pre>def load_sac_eff(...):     # Loads SAC/Efficiency correction in wkspc     "Corr" or sets it to "1"     newCalFile = EXPR_FILE.CorrVanDir + '/' +     EXPR_FILE.GrpFile  def focus_one(...):     newCalFile = join(EXPR_FILE.user,     EXPR_FILE.GrpFile)  def create_vana(EXPR_FILE, NoAbs=False):     if EXPR_FILE.VGrpfocus == "sam":         GrpFile = EXPR_FILE.Path2DatGrpFile     else:         GrpFile = EXPR_FILE.Path2VanGrpFile     . . .     DiffractionFocussing(InputWorkspace="Vanadium_corr",     OutputWorkspace="Vanadium_foc",     GroupingFileName=GrpFile,     PreserveEvents=False)</pre>
VrunnoList / Vanadium	39191	-
VERunnoList / V-Empty	39187	-
SERunnoList / S-Empty	0	-
ExistV / ExistingV	Yes	<pre>def focus_all(...)  if EXPR_FILE.ExistV == "load":     for i in EXPR_FILE.bankList:         spec = i - 1         vanfil = EXPR_FILE.CorrVanFile + "-" +         str(spec) + ".nxs"         LoadNexusProcessed(Filename=vanfil,         OutputWorkspace="Vanadium-" + str(i))         # CORRECT elif EXPR_FILE.ExistV == "no" and EXPR_FILE.VGrpfocus == "van":     print "was here?"     cry_vana.create_vana(EXPR_FILE, NoAbs=NoVabs)</pre>
XYEtof / XYE-TOF	Yes	-
XYEdspc / XYE-D	Yes	-
GSS	yes	<pre>def rearrang4gss(OutputFile, EXPR_FILE):     if EXPR_FILE.GSS == "no":         return     if len(EXPR_FILE.bankList[1:]) &gt; 1:         SaveGSS(InputWorkspace="ResultTOFgrp")     else:         SaveGSS(InputWorkspace="ResultTOF-1")</pre>

Nex	Yes	<pre>def rearrange_4nex(OutputFile, EXPR_FILE):     if EXPR_FILE.Nex == "no":         return     SaveNexusProcessed(Filename=OutputFile + ".nxs", InputWorkspace="ResultTOFgrp")</pre>
NBank	3	<pre>def sets_drangle(wkspc, EXPR_FILE):     datamatrix = mtd[wkspc]     for i in range(0, EXPR_FILE.Nbank):  def __init__(self, instr):     self.nbank = NUM_BANK_DIC[self.sname]  def load(...) blist = range(1, self.Instr.nbank + 1)</pre>
CropRange / CropRange	0.05 0.95	<pre>def sets_drangle(wkspc, EXPR_FILE):     datamatrix = mtd[wkspc]     for i in range(0, EXPR_FILE.Nbank):         x_data = datamatrix.readX(i)         last = len(x_data) - 1         CropRange = EXPR_FILE.CropRange[i].rstrip().split()         xbegin = str(x_data[0] * (1 + float(CropRange[0])))         xend = str(x_data[last] * float(CropRange[1]))         datbin = math.exp(math.log(x_data[last] / x_data[0]) / last) - 1         if datbin &gt; float(EXPR_FILE.Bining[i]):             print 'WARNING: Rebinning in *pref file ' + EXPR_FILE.Bining[             i] + ' is lower than diffraction focusing rebinning step'             print 'WARNING: Rebinning Kept to be ' + str(datbin) + ' for bank ' + str(i + 1)             EXPR_FILE.Bining[i] = str(datbin)             Drange = xbegin + "," + EXPR_FILE.Bining[i] + "," + xend             EXPR_FILE.Drange.append(Drange)             EXPR_FILE.dataRangeSet = True         return sets_drangle</pre>
Bining / Bining	0.0001	
bankList / BankList	1-3	<pre>def focus_one(...) cry_load.bin_bank("ResultD", EXPR_FILE.bankList, EXPR_FILE.Drange) . . . cry_load.bin_bank("ResultD", EXPR_FILE.bankList, EXPR_FILE.Drange)     for i in EXPR_FILE.bankList:         ConvertUnits()         ReplaceSpecialValues()         ReplaceSpecialValues()  def divide_samp_vana(EXPR_FILE, Norm):     for i in EXPR_FILE.bankList:         RebinToWorkspace()         Divide()     else:         for i in EXPR_FILE.bankList:             RenameWorkspace()</pre>

		<pre>def strip_the_vana(EXPR_FILE, LoadUnstrip=""):     cry_load.split_bank("Vanadium", bankList=EXPR_FILE.bankList, Del=True)</pre> <pre>def rearrange_4xye(...):     if (units=="D" and EXPR_FILE.saveXYEd) or (units == "TOF" and EXPR_FILE.saveXYEtof):         for i in EXPR_FILE.bankList:             inwkpsc = "Result" + units + "-" + str(i)             SaveFocusedXYE()</pre> <pre>def split_bank(InputArea, bankList, Del=True):     for i in bankList:         CropWorkspace()     if Del:         mtd.remove(InputArea)</pre> <pre>def bin_bank(InputArea, bankList, Drange):     for i in bankList:         Rebin()</pre> <p>Many other places too... mostly as a loop in which which n algorithm is applied</p>
--	--	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------