Powder Diffraction Reduction Status Update 16.11.17, ILL

News since last meeting

Large refactor of the calibration algorithm

- Better organization of the code
 Important speed-up (~1min max)
 Interpolation of overlapping 2theta angles
 ROI normalisation for calibration
- Fixed the exclusion regions for calibration

Open questions

Rebinning of the temperature axis

- Currently Rebin is used which is designed for histogram data
- It is not correct for non-equidistant points
- How exactly we want to group our points?
- New algorithm GroupPointData?

Picking up the scanning observable

- Currently this has to be typed by the user
- User should choose e.g. from [sample.temperature, Omega.Position]
- Can we define the most common parameters, and suggest a list?
- Can we make this fully automatic?


Detector scan data reduction

- Currently detector scan reduction logic is contained in the calibraiton algorithm
- Should we rename it to PowderDiffDetectorScan, with options [ProduceCalibration, Reduce]?
- And/or make 2 separate light-weight wrappers?


A side note about dead cell definition

- Few cells count 0, but not always
- A cell is treated as dead, if it counts zero in 80% or more cases

Investigated the difference in calibration constants for few cells




The ratio histogram for cell #162



Outside the beam-stop it is indeed about 0.7

It turned out that this comes from raw data.



Some cells do count too much, and get a smaller calibration constant.

23/11/2017, 08:38 reveal.js

Status

- PR is open, ready for review
- New package is built
 Built on the server this time