


# Data reduction workflow

**Note:** Dashed boxes are optional, a dashed path can be overridden by user input

Dataparams



## Data Worksheet

The Data Reader GUI should have the appearance of a spreadsheet (e.g. LAMP), so that the user has at one glance an overview of the required information (sample and background scattering/transmission, dark current) and already given run numbers. The spreadsheet should provide space for complementary instrument configuration (e.g. detector positions) and should allow to add new rows to append several samples.

In general, each run can contain parameter frames (axis/wavelength/time slices). Features like adding series of runs (summing intensities +total acquisition time/monitor of equivalent frames) as well as concatenating runs in parameter frame depth have to be supported.

## Spin-leakage correction

As a prerequisite for spin leakage correction, data (optionally user) has to provide a flag for the experiment [polarised beam only (SANSPOL) or polarisation analysis (POLARIS)] and the order of spin states in parameter depth. Algorithms for spin leakage correction for SANSPOL (i.e. 2 spin states) and POLARIS (4 spin states) can be found in /grasp/polarisation analysis/sanspol\_correct.m and pa\_polarisation\_correct.m respectively. Note that using a 3He cell the correction factors are time-dependent. These factors are given by the time evolution of the transmission without sample in beam.

## Visualisation

Basic functionalities for graphical output should comprise:

Data point selector/cursor

Zoom function

Select displayed detector area as ROI or mask

For processed 2D detector data:

Scroll through depth

Movie function in depth (e.g. along time dimension or instrument parameter)

Optional: store results as annotation in Nomad logfile (logs.ill.fr)

