



POWDER DIFFRACTION REDUCTION


GAGIK, 22.09.2017, ILL



CURRENT STATUS

- ✓ Instrument geometry
- ✓ Loader validated
- ✓ Scanning detector loading
- ✓ Reduction algorithm
- ✓ Calibration algorithm
- ✓ Fullprof export
- ❑ Code cleanup and speedup
- ❑ Unit tests and documentation
- ❑ ROC correction
- ❑ GSAS export
- ❑ Omega scan
- ❑ DIB case

✓ LOADING VALIDATED



After final fix to the $2\theta_0$ rotation treatment and dead panels allocation (both on the right end), obtained a good match with LAMP in signed 2θ axis.

REDUCTION AND CALIBRATION ALGORITHMS

PowderDiffILLReduction input dialog

Performs powder diffraction data reduction for ILL instrument D20.

Run TEST/RC_171/967420.nxs

CalibrationFile /Desktop/D20_median.nxs

ROCFile

NormaliseTo Monitor

ROI 0,153.6

Observable sample.temperature
Scanning observable, a Sample Log entry

☐ SortObservableAxis

Unit ScatteringAngle

OutputWorkspace

? Keep Open ☐

PowderDiffILLCalibration input dialog

Performs detector efficiency correction calculation for powder diffraction instrument D20 at ILL.

CalibrationRun EST/d20_calib/967076.nxs

CalibrationFile

CalibrationMethod Median

NormaliseTo None

ROI 0,153.6

ExcludedRange -3.2,0

PixelRange 65,3072


OutputResponseWorkspace

OutputWorkspace median


? Keep Open ☐

CALIBRATION CONSTANTS

Calibration derived in LAMP




Calibration derived in Mantid, median method



- ✓ Similar results
- Note the different axes
- Note the different methods


CALIBRATION CONSTANTS



- ✓ Mean, median and max likely methods implemented
- ✓ Good (and similar) results
- ✓ Takes about 1m 30s for median and mean
 - ❑ Maximum likelihood too slow (45m) in python
 - ❑ Investigate multiprocessing or C++?
- ✓ Monitor normalisation
 - ❑ Is time normalisation needed?
 - ❑ What about ROI normalisation?
 - ❑ Exclude dead pixels from final absolute scaling?

Currently no interpolation can be done
Assuming the scan step is integer multiple of the pixel size


CALIBRATED DATA COMPARISON




✓ Calibration works

CALIBRATION CLOSURE - DRIFT


Calibration closure - drift of 1.5%



COMBINED RESPONSE FROM VANADIUM




- Optional output of the calibration algorithm as a control plot
- Contains counts of all the pixels iteratively weighted-summed in the overlapping regions



- Steep fall at the verge of the scan
- Present also in raw data
- Not an artifact from calibration procedure

FULLPROF EXPORT OPERATIONAL



- ✓ **SaveFocusedXYE** for FullProf format 10
- ☐ Need to add the temperature in the header of the file
- ☐ Need to crop out the dead pixels

TO DO

- ☐ Wrap up D20 calibration
- ☐ Build a package from the branch for testing
- ☐ Write unit tests, documentation, workflow diagrams, documentation tests and usage examples.
- ☐ Any other feature to include in the algorithms?
- ☐ GSAS export
 - ☐ Need to figure out how to load existing .gsa or .gda to GSAS II, then try to mimic the format
 - ☐ Research if the formats provided by SaveGSS existing in Mantid satisfy our case
 - ☐ Implement new format saving otherwise
- ☐ ROC Correction – **numors needed**
- ☐ Omega scan – **numors needed**
- ☐ DIB – **numors needed**